



## VITREX® Steel GENERAL PURPOSE CAPILLARY BLOOD LANCET

- ◆ LOW COST simple, flat, steel lancet
- ◆ SAMPLE VOLUME excellent blood flow
- ◆ HIGH QUALITY hardened surgical steel
- ◆ MINIMUM DISCOMFORT modern asymmetric design
- ◆ STERILE sterilised by gamma radiation to DIN EN 552/556

## PRODUCT

PRODUCT	DESCRIPTION	ARTICLE NUMBER
VITREX® Steel	Flat blade sterile lancet	350213

## PACKAGING

PRIMARY BOX	CASE	SHIPPER	PALLET
200 pcs	50 x Primaries (10,000pcs)	4 x cases (40,000pcs)	720,000 pcs

## PRODUCT SPECIFICATION

BLADE	
MATERIAL	Stainless Steel Grade X 5 Cr18Ni9 to DIN 17440
TOTAL LENGTH	49.00 mm (+/- 0.50 mm)
EXPOSED LENGTH	3.0 mm (+/- 0.30 mm)
POINT	12° asymmetric point with 18° angle

  

STERILITY	
METHOD	Gamma irradiation by cobalt 60 source to DIN EN 552/556
SHELF LIFE	5 years
WRAPPING PAPER	
TOP (PRINTED)	Supraferm 40g/m <sup>2</sup> kraft paper 20g/m <sup>2</sup> polyethylene coated
BOTTOM (TRANSPARENT)	Supraferm 30g/m <sup>2</sup> bleached pergamyne 20g/m <sup>2</sup> polyethylene coated